











# INSTRUCTION FOR USE

## SwabSURE ListeriaP

Samples should be kept cool until testing is performed.

<p><b>STEP 1</b></p> <p>Identify the site to be sampled. It is recommended to use a template for accurate results where possible (TS/15-T).</p> 	<p><b>STEP 2</b></p> <p>Hold the swab with both hands, pull down and twist the swab cap away from the tube; breaking the security seal in the process.</p> 
<p><b>STEP 3</b></p> <p>Remove swab from tube using the cap as a holder. Do not touch the swab tip or shaft and make sure that it does not touch any surface before sampling. <b>Do not discard the tube.</b></p> 	<p><b>STEP 4</b></p> <p>Sample designated area using a template where possible. If the area permits, swab in 3 planes (horizontal, vertical and diagonal) whilst rotating swab between forefinger and thumb.</p> 
<p><b>STEP 5</b></p> <p>When sampling is completed insert swab back into the tube and push cap on tightly.</p> 	<p><b>STEP 6</b></p> <p>Write sample details on label. Store at 4 - 8°C (or cool &amp; dry as possible) and send to the laboratory within 24 hours for chromogenic test.</p> 
<p><b>Chromogenic Test</b></p> 	<p><b>STEP 7</b></p> <p>Alternatively (where permitted), immediately insert swab into SS-L01 media tube. Pull down on the shaft of the swab whilst supported inside the tube to activate breakpoint.</p> 
<p><b>STEP 8</b></p> <p>Place tubes into a rack &amp; incubate them horizontally at 37°C for 24-48 hours</p> 	<p><b>STEP 9</b></p> <p>Check results in 24-48 hours. Change of media colour to turquoise-blue indicates presence of <i>L.monocytogenes</i> or <i>L.ivanovii</i></p> 

Have you considered? Sterile Sampling Template - TS/15-T